

Bedding in a Box

HCPCS:



FIGURE A

- The Bedding in a Box set is perfect for use with most standard manual, semi-electric and full electric beds
- Set contains 1 pillow case, 1 fitted sheet, 1 cover sheet and 1 blanket that are ideal for standard and hospital homecare beds Fitted sheet measures 36" x 80" x 5" - 6"
- Sheets are made of a high quality cotton and polyester blend Blanket is 100% cotton
- All components are machine washable

PRODUCTS	Description
15030HBC	Bedding in a Box,1 EA,5/CS