

## Hospital Bed Linen Kit

HCPCS:



FIGURE A

- The Fitted Sheet is perfect for use with most standard manual, semi-electric and full electric beds
- Contains 2 fitted sheets that are ideal for standard hospital and homecare beds Fitted sheet measures 36" x 80" x 5" - 6"
- Sheets are made of a high quality cotton and polyester blend
- All components are machine washable

### PRODUCTS

### Description

|          |                                   |
|----------|-----------------------------------|
| 15030HBL | Hospital Bed Linen Kit,1 EA,10/CS |
|----------|-----------------------------------|