

Half Length Bed Rail Bariatric Lightweight Beds

HCPCS:



- Bed Rails may be used as an assistive device or for fall prevention for individuals utilizing a medical homecare bed. Half rails 15021BV-LW are designed to be used with Drive's Lightweight Bariatric homecare bed series (Models # 15300LW, 15302LW, and 15303LW). Use of these rails with any other bed system is not recommended. With a limited lifetime warranty, the bariatric bed rails ensure patients are always supported.
- Durable steel construction in rich brown finish offers attractive and easy-to-maintain homecare style
- Easy to install half Rail designed specifically for Lightweight Bariatric Homecare Bed series
- Product Safety & Performance Standards: Incompliance with FDA Hospital Bed System Dimensional Assessment Guidance to Reduce Entrapment as well as CPSC ASTM standards
- Ideal for maximum protection and support when bed is elevated

PRODUCTS	Description
15021BV-LW	Half Side Rail,42in/48/54 LW Baribed,1pr

SPECIFICATIONS	Description
Warranty	Steel Frame: Limited Lifetime Knobs: 1 Year
Color	Brown