

## M-Rail Home Bed Assist Handle

HCPCS:



FIGURE A




FIGURE B

- Designed to provide support when getting into and out of bed and when moving from a lying position to a sitting position
- M shape allows for multiple hand positions to accommodate any user
- Multiple crossbars creates easy grip for any height use
- Easy, quick assembly
- Safety strap wraps around mattress or box spring to ensure a secure fit
- Special extra-long safety strap enables user to be in a single, twin, double, queen or king size bed
- Made of powder coated steel for long life and easy cleaning
- Includes accessory pouch for storing remote controls, glasses, magazines and more

### PRODUCTS

### Description

1222P	M Rail with pouch
-------	-------------------

### SPECIFICATIONS

### Description

Warranty	Limited Lifetime
Weight Capacity	300 lbs