

Multi-Height Manual Bed

HCPCS:




FIGURE A




FIGURE B



FIGURE C



FIGURE D

- Convenient fold down hand cranks at the foot end for Foot, Head and overall bed height adjustment
- Easy to set up, bed ships in two cartons, shaft stores under frame when unassembled
- Head/foot adjustment provide anatomically correct sleep surface
- Reinforced channel frame construction offers superior strength, resists twisting and reduced weight
- Long lasting Zinc-coated spring deck for Pt comfort
- Heavy Duty bed ends with durable HD Polyethylene Molded Bed End covers are attractive, easy-to-maintain
- Satisfies applicable FDA entrapment guidelines
- Patient weight Capacity 350 lbs. (158kg). Total Weight Capacity 450 lbs. (204kg)
- Half rails and full rails available in packages (PKG) or sold separately

PRODUCTS

Description

| | |
|-----------------|-----------------------------|
| 15003P | Manual Bed,1/cv |
| 15003P-FR | Manual Bed w/Full Rails Pkg |
| 15003P-HRPKG-TS | Manual Bed HRail TSMatt Pkg |

SPECIFICATIONS

Description

| | |
|--------------------|--|
| Warranty | (Original Purchaser) Limited Lifetime on Welds, 5 Years on Frame, 1 Year on All Other Parts and Components |
| Suitable Room Size | 36" (.914m) (W) x 88" (2.24m) (L) |
| Weight | 156 lbs. (71kg) |