

## Semi-Electric Bed (Single Crank)

HCPCS:




FIGURE A



FIGURE B

- Quiet, smooth operation, Head/Foot section adjusts electronically, Height adjusts manually
- Easy-to-transport and set up, Bed ships in two cartons
- Shaft stores under frame when unassembled
- Channel frame construction provides superior strength and reduced weight
- Reinforced frame resists twisting and bending.
- Long Lasting Zinc-coated spring deck for Pt comfort
- Heavy Duty bed ends with durable HD Polyethylene Molded Bed End covers are attractive, easy-to-maintain
- Clip-on pendant mounts on bed rail, operate Head and Foot section simultaneously
- Half rails and full rails available in packages (PKG) or sold separately
- Satisfies applicable FDA entrapment guidelines
- Extension kit available (model #15005EXTKIT-L) extends length from 80" to 84"
- Patient weight Capacity 350 lbs. (158kg). Total Weight Capacity 450 lbs. (204kg)
- Motor is UL Approved, Bed is manufactured in ISO9001 & ISO13485 facility

| PRODUCTS        | Description                   |
|-----------------|-------------------------------|
| 15004P          | SemiElec Bed,1/cv             |
| 15004P-FR       | SemiElec Bed w/Full Rails Pkg |
| 15004P-FRPKG    | SemiElec Bed,FRail IMatt Pkg  |
| 15004P-HR       | SemiElec Bed w/Half Rails Pkg |
| 15004P-HRPKG    | SemiElec Bed HRail IMatt Pkg  |
| 15004P-HRPKG-TS | SemiElec Bed HRail TSMatt Pkg |

| SPECIFICATIONS         | Description                                                                                                               |
|------------------------|---------------------------------------------------------------------------------------------------------------------------|
| Warranty               | (Original Purchaser) Limited Lifetime on Welds, 5 Years on Frame, motor, pendant, 1 Year on All Other Parts and Component |
| Carton Shipping Weight | 176 lbs. (80kg)                                                                                                           |
| Weight                 | 161 lbs. (73kg)                                                                                                           |