

## Semi-Electric Bed (Single Crank)

HCPCS:




FIGURE A



FIGURE B

- Quiet, smooth operation, Head/Foot section adjusts electronically, Height adjusts manually
- Easy-to-transport and set up, Bed ships in two cartons
- Shaft stores under frame when unassembled
- Channel frame construction provides superior strength and reduced weight
- Reinforced frame resists twisting and bending.
- Long Lasting Zinc-coated spring deck for Pt comfort
- Heavy Duty bed ends with durable HD Polyethylene Molded Bed End covers are attractive, easy-to-maintain
- Clip-on pendant mounts on bed rail, operate Head and Foot section simultaneously
- Half rails and full rails available in packages (PKG) or sold separately
- Satisfies applicable FDA entrapment guidelines
- Extension kit available (model #15005EXTKIT-L) extends length from 80" to 84"
- Patient weight Capacity 350 lbs. (158kg). Total Weight Capacity 450 lbs. (204kg)
- Motor is UL Approved, Bed is manufactured in ISO9001 & ISO13485 facility

PRODUCTS	Description
15004P	SemiElec Bed,1/cv
15004P-FR	SemiElec Bed w/Full Rails Pkg
15004P-FRPKG	SemiElec Bed,FRail IMatt Pkg
15004P-HR	SemiElec Bed w/Half Rails Pkg
15004P-HRPKG	SemiElec Bed HRail IMatt Pkg
15004P-HRPKG-TS	SemiElec Bed HRail TSMatt Pkg

SPECIFICATIONS	Description
Warranty	(Original Purchaser) Limited Lifetime on Welds, 5 Years on Frame, motor, pendant, 1 Year on All Other Parts and Component
Carton Shipping Weight	176 lbs. (80kg)
Weight	161 lbs. (73kg)