

## Prime Care Bed Model P903

HCPCS:



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E




FIGURE F




FIGURE G




FIGURE H

- 600 lbs. weight capacity\*
- 7" super low height can help minimize falls (Figure A) while 30" high working height can facilitate staff &
- Integrated width with push button expansion from 36"-39"-42" provides extra measure of comfort and dignity (Figure I)
- 80", 84" and 88" sleep surface length helps comfortably accommodate the tallest residents
- 10 function Backlit Hand Pendant includes Auto Contour, Cardiac Comfort Chair (Figure D) and under bed night light provides an added measure of security at night
- Synchronized Hi/Lo control system keeps bed level for resident comfort and positioning
- Welded Steel construction, slat support deck is durable and easy to clean
- Roll at any height convenience with 3" heavy-duty dual footprint casters plus steer control facilitate easy maneuverability (Figure H)
- Color coded, dual pedal floor lock system keeps bed securely in place when engaged (Figure G)
- Adjustable wall bumper helps protect walls and bed from damage (Figure E)
- Standard features include: adj. mattress retainers, heel lift, IV pole holders
- Optional: bed ends (two styles and three wood grain laminate colors), Assist bars, rails, deluxe assist rails, trapeze, embedded staff control, battery backup, transporter
- Trendelenburg/Reverse Trendelenburg capable with staff control option (Figure C)
- Warranty: 15 yr. on frame, lifetime on welds, 5 yr. on electronics

| PRODUCTS                 | Description                                 |
|--------------------------|---|
| P903                     | Prime Care Bed Model 903                    |
| SPECIFICATIONS           | Description                                 |
| Warranty                 | Lifetime on Welds                           |
| Wheels                   | 550 lbs                                     |
| Optional Width Expansion | 36"- 39"-42"                                |
| Sleep Surface            | 80", 84", 88"                               |
| Warranty Limited         | 15 Year on Frame                            |
| Warranty on Batteries    | 5 years on Electrical                       |
| ACCESSORIES              | Description                                 |
| PLTCAR                   | LTC Assist Rail, 1PR, PR                    |
| PLTCAB                   | LTC Assist Bar, 1PR, PR                     |
| PLTCDR                   | LTC Deluxe Rail, 1PR, PR                    |
| P903SC05                 | P903 Staff control for Ti-Motion Motors, EA |
| P903SC03                 | P903 Staff control for Linak Motors, EA     |
| 42BWBE-A                 | 42" BRNWOOD BEDEND, CHERRY, EA, EA          |

## ACCESSORIES

## Description

|            |  |
|------------|--|
| 42MCBE-FM  | 42" MILL CRK BED END,FIG MAHOG,EA, EA        |
| 35MCBE-A   | 35" MILL CRK BEDEND,CHERRY,EA, EA            |
| 42BWBE-B   | 42" BRNWOOD BEDEND,OKOAK,EA, EA              |
| 35BWBE-A   | 35" BRNWOOD BEDEND,CHERRY,EA, EA             |
| 35MCBE-B   | 35" MILL CRK BEDEND,OKOAK,EA, EA             |
| 35BWBE-B   | 35" BRNWOOD BEDEND,OKOAK,EA, EA              |
| P903-T-BBK | BATTERY BACKUP W/MOUNTING, P903, Timotio, EA |