

## Primus LTC Bed Scale with Hardware

HCPCS:



- 1,000-pound weight capacity
- LED display for easily readable measurements
- Low battery indicator
- "Weigh/Zero" button to zero scale before use
- Quickly change between pound and kilogram with one touch of a button
- Clearly labeled head and foot sections
- "MEM/Recall" button freezes displayed weight and stores in memory
- Lightweight and portable to connect to Drive's P750 Bed

### PRODUCTS

### Description

PBEDSCALE-PKG

Primus LTC Bed Scale w/ hrdwr