

C/D/E Cylinder Rack, 6 Cylinder Capacity

HCPCS:

- The C/D/E Cylinder Rack ensures durability and strength in storing your oxygen cylinders with an all-welded frame construction.
- It also features individual spacing, positive cylinder placement, and weight distribution to promote safe storage.
- Features an all-welded frame construction to ensure durability and strength
- Promotes safe storage with individual spacing, positive cylinder placement, and weight distribution
- Equipped with bolt down bars and feet
- Chrome finish
- Fits up to 6 C, D, or E cylinders
- Sold individually
- To safely store your oxygen cylinders, purchase the C/D/E Cylinder Rack.

PRODUCTS	Description
CR-EPC6	D/E OXYGEN CYLINDER RACK-6 CYL