

Cylinder Holder for Wheelchair, D/E Cylinders

HCPCS:

- The Cylinder Wheelchair Holder provides a convenient and stable way to store or transport your oxygen cylinder.
- It mounts easily to the back of most wheelchairs with easy-to-assemble U-bolts and holds one D or E cylinder.
- Constructed of durable and sturdy steel
- Versatile 2-in-1 design allows for use on most wheelchairs or as a wall mount
- Bracket easily mounts to back of wheelchair with easy-to-assemble U-bolts
- Designed to hold one (1) D or E oxygen cylinder
- Shown on wheelchair with E cylinder; cylinder NOT included
- Quantity: 4/case

PRODUCTS	Description
150-0110	WHEELCHAIR O2 E TANK HOLDER 4EA/CV