

Drop Arm Steel Commode

HCPCS: E0165,E0165



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E



FIGURE F



FIGURE G



FIGURE H

Drop arm for ease of transfers comes standard with padded seat and deepest seat

PRODUCTS

PRODUCTS	Description
11125PSKD-1	Commode,DropArm,Padded Seat,KD,1/cs

Description

SPECIFICATIONS

SPECIFICATIONS	Description
Product Weight	21 lb
Seat Dimensions	15.25" x 16.75"
HCPCS	E0165
Height Adjustment	19" - 23"
Width Between Arms	19"
CA Prop 65	no