

## Folding Steel Commode Combo with Absorbent Commode Liners

HCPCS:



FIGURE A



FIGURE B



FIGURE C




FIGURE D




FIGURE E



FIGURE F



FIGURE G



FIGURE H

- Includes RTL12095 Super Absorbent Commode Liner 7 Pack
- Easily opens and folds
- Folds flat for convenient storage and transportation
- Tool-free set-up
- Durable plastic snap on seat and lid
- Comes complete with 7.5 qt commode bucket, cover and splash shield
- Grey powder coated steel welded construction increases strength and durability
- Plastic arms for added comfort
- Limited Lifetime Warranty
- Piggy back label supplied on product hangtag
- Made using recycled materials

PRODUCTS	Description
11148-PKG-1	Folding Commode Liner Combo, 1cs
SPECIFICATIONS	Description
Warranty	Limited Lifetime
Seat Dimensions	13.5" (W) x 15" (D)
Weight Capacity	350 lbs
Carton Shipping Weight	13.5 lbs
Weight	13.5 lbs
CA Prop 65	[WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to <a href="http://www.P65Warnings.ca.gov">www.P65Warnings.ca.gov</a> .]