

EasyPulse 5 Conserving Device

HCPCS:

- The EasyPulse 5 is the smallest and lightest oxygen conserving device, weighing in under 10 oz.
- Single dial control and ergonomically designed T-Handle/Hand Tight Nut make setting and connecting the EasyPulse 5 easy for any patient or caregiver.
- With a conserving ratio of up to 5.7:1, the EasyPulse 5 reduces deliveries of cylinders to patients and its proven pneumatic reliability, robust design and lack of batteries helps to reduce overall operating costs.
- Inlet Pressure Range: 300-3000 psig
- Pulse Settings: 1, 2, 3, 4 and 5
- Flow Settings: 2 LPM Continuous
- Cannula Requirements: Single Lumen with up to 7' Tubing
- Operating Altitude Range: Sea Level to 10,000'

PRODUCTS	Description
198705	PRECISION EASYPULSE 5 CONSERVER