

Precision EasyPulse 5 + 6 Oxygen Conserver

HCPCS:

- The Precision EasyPulse 5 + 6 Oxygen Conserving Regulator conserves oxygen and regulates pressure for patients receiving oxygen therapy with a high-pressure cylinder.
- It offers 6 continuous flow settings and 5 pulse settings.
- 6 continuous flow settings of 1, 2, 3, 4, 5, and 6 liters per minute
- Delivers pulse at settings 1-5 like an electronic conserver without the need for batteries
- Delivers a front-loaded tapered wave form; all oxygen is delivered in the first third of the inspiratory phase
- Lightweight at only 9.9 ounces and compact at only 4.75"
- Conserving ratio of up to 5.7 to 1
- Fits all CG-870 post oxygen tanks
- For a versatile, lightweight, and effective oxygen conserving device, purchase the Precision EasyPulse 5+6 Oxygen Conserving Regulator.

PRODUCTS	Description
198712	PRECISION EASPULSE 5 + 6 CONSERVER