

UltraTENS II Leadwire - Single 1EA/

HCPCS:

PRODUCTS	Description
77718	CHATTANOOGA INTELECT IF 1EA/
ACRL-9100CA	ACCURELIEF WIRELESS 4EA/CV
CC3001	Black TENS/EMS Carrying Case 32EA/CV
DI0001X	Intensity Series AC Adapter, 1 1EA/
DI2000-AC	Charger and USB for 2nd Gen. 1EA/
DT2020-BC	Venti TENS Replacement Battery 1EA/
DT7202BC	Replacement Batter Cover for 1EA/
DT7202-ELEC204032C	TENS 7000 2" x 2" 56EA/CV
DU1025X	AC Adapter for US 1000 3rd ed 1EA/
DU2035X	AC Adapter for US Pro 2000 1st 1EA/
LC2828	TheraSonic Conductive Gel, 2-o 1EA/
LL2003	TENS Lotion - 2oz 140EA/CV
LS2177-2	2 oz Conductive Spray 128EA/CV
LS5255CA	TheraSonic Conductive Gel, 5L 1EA/
LS5266CA	TheraSonic Conductive Gel, 16- 12EA/CV
LS5288CA	TheraSonic Conductive Gel, 8-o 24EA/CV
TA2000	Lithium Ion Battery - 2nd Gen 1EA/
TA2020	TENS Cell - AA Batteries 1EA/
WS3147	Snaps-Black with Pigtail Qty: 1EA/
WW3005	Lead Wire - Standard 45 1EA/
WW3766	Lead Wire Extender 1EA/
WW5000	5 Foot Lead Wire(2 per pkg) 1EA/