

AC Adapter for US 1000 3rd Edition

HCPCS:

- The US Pro 1000 A/C Adapter is the replacement power cord for all 3rd Edition US Pro 1000 ultrasound devices.
- It is lightweight and portable, weighing less than a pound.
- Replacement A/C power adapter reliably powers the portable US Pro 1000 ultrasound devices
- Lightweight and portable to be taken everywhere with portable device
- Input: AC 100-240V, 50-60Hz
- Output: 12V 1A 12W
- Compatible with 3rd edition US Pro 1000 devices only
- Sold individually

PRODUCTS

Description

DU1025X	AC Adapter for US 1000 3rd ed 1EA/
---------	------------------------------------
