

AC Adapter for US Pro 2000 1st and 2nd Edition

HCPCS:

- The US Pro 2000 A/C Adapter is the replacement power cord for all 1st and 2nd Edition US Pro 2000 ultrasound devices.
- It is lightweight and portable, weighing less than a pound.
- NOT for use with 3rd Edition US Pro 2000 devices (#DU3035).
- Replacement A/C power adapter reliably powers the portable US Pro 2000 ultrasound devices
- Lightweight and portable
- Compatible with 1st and 2nd edition US Pro 2000 devices only; not suitable for use with 3rd edition devices
- Sold individually

PRODUCTS

Description

DU2035X	AC Adapter for US Pro 2000 1st 1EA/
---------	-------------------------------------