

## TheraTouch EX4 E-Stim Machine

HCPCS:

- The TheraTouch EX4 is an advanced four-channel electrotherapy system offering a wide range of treatment options in one user-friendly, ergonomic design.
- Equipped with a bilingual menu-driven touchscreen interface, the EX4 guides clinicians step-by-step through treatment protocols to ensure the ideal therapeutic treatment including electrode placement for better treatment outcomes.
- 4 channels of stimulation with the following waveforms:
  - Interferential
  - Premodulated
  - Biphasic
  - Asymmetrical
  - Symmetrical
  - RAAS
  - NMS
  - Russian

### PRODUCTS

### Description

|        |                                    |
|--------|------------------------------------|
| DQ7001 | THERATOUCH EX4 w/Therapy Cart 1EA/ |
|--------|------------------------------------|