

## Digital TENS AA Unit

HCPCS:

- This Digital TENS AA Unit is compact and portable, making it ideal for pain relief and treatments anywhere and everywhere.
- It increases circulation of the blood in the muscle tissue by sending pulses directly into the muscles to relieve pain and treat muscle conditions.
- NOTE: This is a Class II Medical Device and federal law restricts this device to sale by or on the order of a practitioner.
- Rx Required.
- Variable wavelengths, pulse rates and treatment modes and a timer for personalized treatments
- Large print LCD display makes unit easy to use and adjust
- Convenient last-used setting recall
- Comes with color coded lead wires, AA batteries, 4 reusable electrodes, and a carrying case
- 2-year warranty
- For convenient and effective electrode therapy, purchase your Digital TENS AA Unit.

### PRODUCTS

### Description

|        |                                     |
|--------|-------------------------------------|
| DA1812 | TENS AA - Digital TENS Unit 40EA/CV |
|--------|-------------------------------------|