

Relief Wrap Conductive Brace with 2" x 2" Electrodes, Extra Large

HCPCS:

- The Relief Wrap Conductive Brace helps users receive electrode therapy in hard to reach areas without requiring the help of others.
- It can be customized to fit patient's needs based on treatment protocol, accommodating either 4 electrodes for TENS and EMS therapy or 8 electrodes for interferential therapy.
- Comes with 4-pack of attachable 2" x 2" electrodes (#EP2020RW2)
- Size extra-large, fits approximately 40" to 45" waists; also available in small, medium, and large
- Easy to maintain, hand washable
- Compatible with most devices
- To independently receive electrode therapy in even the hardest to reach areas, purchase the Relief Wrap Conductive Brace.

PRODUCTS	Description
GF5130	Relief Wrap - XLarge 2X2 1EA/