

Relief Wrap Conductive Brace with 2" x 2" Electrodes, Small

HCPCS:

- The Relief Wrap Conductive Brace helps users receive electrode therapy in hard to reach areas without requiring the help of others.
- It can be customized to fit patient's needs based on treatment protocol, accommodating either 4 electrodes for TENS and EMS therapy or 8 electrodes for interferential therapy.
- Comes with 4-pack of attachable 2" x 2" electrodes (#EP2020RW2)
- Size small, fits approximately 20" to 25" waists; also available in medium, large, and extra large
- Easy to maintain, hand washable
- Compatible with most devices

PRODUCTS	Description
GF5350	Relief Wrap - Small 2x2 1EA/