

Economy Self-Adhesive Electrodes, 2"x 2" Tan Cloth in Poly Bag

HCPCS:

- The Economy Tan Cloth 2" x 2" Electrodes are made with the highest clinical standards and are for use with TENS and EMS units.
- These highly conductive electrodes deliver more current to better distribute across entire pad.
- They are self-adhesive and easy to use, and with premium tan cloth and a full-coverage wire connector, they are durable and safe for any user.
- Self-adhesive pad for ease of use
- Flexible carbon layer allows for maximum charge delivery and comfort on a variety of application sites
- Durable and comfortable premium tan cloth backing
- Durable and reliable full-coverage wire connector
- 4 electrodes/pack; comes sealed in tamper-proof poly bag
- Can be reused up to 15 times per pack
- Purchase the Economy Tan Cloth 2" x 2" Electrodes for effective yet economic electrotherapy.

PRODUCTS	Description
E1P2020TC2	Elec Poly 2x2 50EA/CV/400CT