

Rubber Electrodes, 2.5" x 3.5", Polybag

HCPCS:

- The 2.5" x 3.5" Professional Rubber Electrodes are made with the highest clinical standards and are for use with TENS and EMS units.
- They are made with premium rubber and a full-coverage wire connector, they are durable and safe for any user.
- Flexible and conformable to fit comfortably on a variety of application sites
- Durable, latex-free rubber backing
- Durable and reliable full-coverage wire connector; fits standard 2mm lead wires
- 4 electrodes/pack; comes sealed in tamper-proof polybag
- Size: 2.5" x 3.5"
- Can be reused up to 15 times per pack

PRODUCTS	Description
ER2535B2	Rubber Electrodes, 2.5 x 3.5 1EA/