

Self-Adhesive Electrodes, 1.5" x 14" White Foam, Poly Bag

HCPCS:

- The Rectangular White Foam Electrodes are made with the highest clinical standards and fit most TENS, EMS, IF, and other electro-stim applications.
- These highly conductive electrodes deliver more current to better distribute across entire pad.
- With premium white clinical foam and two full-coverage wire connectors, they are durable and safe for any user.
- Self-adhesive pad for ease of use
- Flexible carbon layer allows for maximum charge delivery and comfort
- Durable and comfortable white foam backing
- Two (2) durable and reliable full-coverage wire connectors
- Size: 1.5" x 14"
- Quantity: 1 electrode/pack
- Can be reused up to 15 times per pack
- For effective and high-quality electrotherapy, purchase the Rectangular White Foam Electrodes.

PRODUCTS

Description

EP15140WF2	Elec Poly 1.5x14 200EA/CV/800CT
------------	---------------------------------