

Self-Adhesive Electrodes, 2" Round White Foam, Foil Pouch

HCPCS:

- The 2" Round White Foam Electrodes are made with the highest clinical standards and are for use with TENS and EMS units.
- These highly conductive electrodes deliver more current to better distribute across entire pad.
- They are self-adhesive and easy to use, and with premium white foam and a full-coverage wire connector, they are durable and safe for any user.
- Self-adhesive pad for ease of use
- Flexible carbon layer allows for maximum charge delivery and comfort
- Durable and comfortable premium white foam backing
- Durable and reliable full-coverage wire connector
- Tyco gel is self-adhesive and made in the USA
- 4 electrodes/pack; comes sealed in tamper-proof foil bag
- 10 packages/foil bag
- Size: 2" round

PRODUCTS

Description

EF2000WF2	Elec Foil 2 Round WF 50EA/CV/400CT
-----------	------------------------------------