

## Self-Adhesive Electrodes, 2" x 2" Blue Cloth in Poly Bag

HCPCS:

- The Blue Cloth Electrodes are made with the highest clinical standards and fit most TENS, EMS, IF, and other electro-stim applications.
- These highly conductive electrodes deliver more current to better distribute across entire pad.
- They are self-adhesive and easy to use, and with premium blue cloth and a full-coverage wire connector, they are durable and safe for any user.
- Self-adhesive pad for ease of use
- Flexible carbon layer allows for maximum charge delivery and comfort on a variety of application sites
- Durable and comfortable premium blue cloth backing
- Durable and reliable full-coverage wire connector
- 4 electrodes/pack
- Size: 2" x 2"
- Can be reused up to 15 times per pack
- Purchase Blue Cloth Electrodes for effective yet economic electrotherapy.

<b>PRODUCTS</b>	<b>Description</b>
E1P2020BC2	Electrode - Poly - 2x2 400EA/CV