

Self-Adhesive Electrodes, 2" x 2" White Foam with Tyco Gel in Foil Pouch

HCPCS:

- The White Foam Electrodes are made with the highest clinical standards and fit most TENS, EMS, IF, and other electro-stim applications.
- These highly conductive electrodes deliver more current to better distribute across entire pad.
- With premium white clinical foam and a full-coverage wire connector, they are durable and safe for any user.
- Self-adhesive pad for ease of use
- Flexible carbon layer allows for maximum charge delivery and comfort on a variety of application sites
- Durable and comfortable white foam backing
- Durable and reliable full-coverage wire connector
- Tyco gel is self-adhesive and made in the USA
- Sealed in a tamper-proof foil pouch
- Size: 2" x 2"
- Quantity: 4 electrodes/pack, 10 packages/bag

PRODUCTS

Description

EF2020WF2	Elec Foil 2x2WF 50EA/CV/400CT
-----------	-------------------------------