

Self-Adhesive Electrodes, 2" x 3.5" Tan Cloth with Tyco Gel, Foil Bag

HCPCS:

- The Tan Cloth Electrodes are made with the highest clinical standards and fit most TENS, EMS, IF, and other electro-stim applications.
- These highly conductive electrodes deliver more current to better distribute across entire pad.
- With premium tan cloth and a full-coverage wire connector, they are durable and safe for any user.
- Self-adhesive pad for ease of use
- Flexible carbon layer allows for maximum charge delivery and comfort on a variety of application sites
- Durable and comfortable premium tan cloth backing
- Durable and reliable full-coverage wire connector
- Tyco gel is self-adhesive and made in the USA
- Size: 2" x 3.5"
- Sealed in a tamper proof foil bag
- 4 electrodes/pack, 10 packs/bag

PRODUCTS

Description

EF2035TC2

Electrode - Foil - 2x3.5 - Tan 400EA/CV
