

## Self-Adhesive Electrodes, 2" x 3.5" White Cloth

HCPCS:

- The White Cloth Electrodes are made with the highest clinical standards and fit most TENS, EMS, IF, and other electro-stim applications.
- These highly conductive electrodes deliver more current to better distribute across entire pad.
- With premium white cloth and a full-coverage wire connector, they are durable and safe for any user.
- Self-adhesive pad for ease of use
- Flexible carbon layer allows for maximum charge delivery and comfort on a variety of application sites
- Durable and comfortable white cloth backing
- Durable and reliable full-coverage wire connector
- Tyco gel is self-adhesive and made in the USA
- Sealed in a tamper-proof foil pouch
- Size: 2" x 3.5"
- Quantity: 4 electrodes/pack, 10 packages/bag

<b>PRODUCTS</b>	<b>Description</b>
EF2035WC2	Electrode - Foil - 2x3.5 - 400EA/CV