

## Self-Adhesive Electrodes, 2" x 4" Oval Tan Cloth, Foil Pouch

HCPCS:

- The 2" x 4" Oval Tan Cloth Electrodes are made with the highest clinical standards and are for use with TENS and EMS units.
- These highly conductive electrodes deliver more current to better distribute across entire pad.
- They are self-adhesive and easy to use, and with premium tan cloth and a full-coverage wire connector, they are durable and safe for any user.
- Self-adhesive pad for ease of use
- Flexible carbon layer allows for maximum charge delivery and comfort on a variety of application sites
- Durable and comfortable premium tan cloth backing
- Durable and reliable full-coverage wire connector
- Tyco gel is self-adhesive and made in the USA
- 320/CTN, each sealed in a tamper-proof foil bag
- Size: 2" x 4" Oval
- Can be reused up to 15 times per pack

### PRODUCTS

### Description

---

|           |                                       |
|-----------|---------------------------------------|
| EF2040TC2 | Electrode - Foil - 2x4 - Tan 400EA/CV |
|-----------|---------------------------------------|

---