

Self-Adhesive Electrodes, 2" x 4" White Foam Oval, Foil Pouch

HCPCS:

- The 2" x 4" White Foam Oval Electrodes are made with the highest clinical standards and are for use with TENS and EMS units.
- These highly conductive electrodes deliver more current to better distribute across entire pad.
- They are self-adhesive and easy to use, and with premium white foam and a full-coverage wire connector, they are durable and safe for any user.
- Self-adhesive pad for ease of use
- Flexible carbon layer allows for maximum charge delivery and comfort
- Durable and comfortable premium white foam backing
- Durable and reliable full-coverage wire connector
- Tyco gel is self-adhesive and made in the USA
- Size: 2" x 4" oval
- Quantity: 400/ctn, each sealed in a tamper-proof foil pouch
- Can be reused up to 15 times per pack

PRODUCTS	Description
EF2040WF2	Elec Foil 2x4 WF 40EA/CV/320CT