

## Self-Adhesive Electrodes, 3" Round Tan Cloth, Foil Pouch

HCPCS:

- The 3" Round Tan Cloth Electrodes are made with the highest clinical standards and are for use with TENS and EMS units.
- These highly conductive electrodes deliver more current to better distribute across entire pad.
- They are self-adhesive and easy to use, and with premium tan cloth and a full-coverage wire connector, they are durable and safe for any user.
- Self-adhesive pad for ease of use
- Flexible carbon layer allows for maximum charge delivery and comfort
- Durable and comfortable premium tan cloth backing
- Durable and reliable full-coverage wire connector
- 4 electrodes/pack; comes sealed in tamper-proof foil bag
- Size: 3" round
- Can be reused up to 15 times per pack

### PRODUCTS

### Description

|           |                                    |
|-----------|------------------------------------|
| EF3000TC2 | Elec Foil 3 Round TC 40EA/CV/320CT |
|-----------|------------------------------------|