

## Self-Adhesive Electrodes, 3" Round White Cloth, Foil Pouch

HCPCS:

- The 3" Round White Cloth Electrodes are made with the highest clinical standards and are for use with TENS and EMS units.
- These highly conductive electrodes deliver more current to better distribute across entire pad.
- They are self-adhesive and easy to use, and with premium white cloth and a full-coverage wire connector, they are durable and safe for any user.
- Self-adhesive pad for ease of use
- Flexible carbon layer allows for maximum charge delivery and comfort on a variety of application sites
- Durable and comfortable premium white cloth backing
- Durable and reliable full-coverage wire connector
- Tyco gel is self-adhesive and made in the USA
- 4 electrodes/pack; comes sealed in tamper-proof foil bag
- 10 packages/bag
- Size: 3" round

### PRODUCTS

### Description

EF3000WC2

Elec Foil 3 Round WC 40EA/CV/320CT