

Self-Adhesive Electrodes, 3" Round White Cloth, Poly Bag

HCPCS:

- The 3" Round White Cloth Electrodes are made with the highest clinical standards and are for use with TENS and EMS units.
- These highly conductive electrodes deliver more current to better distribute across entire pad.
- They are self-adhesive and easy to use, and with premium white cloth and a full-coverage wire connector, they are durable and safe for any user.
- Self-adhesive pad for ease of use
- Flexible carbon layer allows for maximum charge delivery and comfort
- Durable and comfortable premium white cloth backing
- Durable and reliable full-coverage wire connector
- 4 electrodes/pack; comes sealed in tamper-proof poly bag
- Size: 3" round
- Can be reused up to 15 times per pack
- Purchase the 3" Round White Cloth Electrodes for effective yet economic electrotherapy.

PRODUCTS

Description

E1P3000WC2	Elec Poly 3 Round White 50EA/CV/400CT
------------	---------------------------------------