

## Parker Labs Signa Spray Electrode Solution, 2-oz

HCPCS:

- The Parker Labs Signa Spray Electrode Solution preps the skin prior to attaching electrodes to maximize their effectiveness and prolong electrode life.
- It is a non-gritty solution that is ideal for use with muscle stimulation equipment and conductive garments, but not recommended for defibrillation.
- Preps the skin to maximize conductivity of electrodes
- Rehydrates self-adhesive electrodes and prolongs electrode life
- Easy to apply and clean; non-gritty with no residue build-up
- Bacteriostatic and hypoallergenic; will not irritate skin
- Will not stain clothes or equipment
- 2-oz size; also available in 8.5-oz size
- Not recommended for defibrillation
- Proudly made in the USA

### PRODUCTS

### Description

---

LS2177-2

2 oz Conductive Spray 128EA/CV

---