

Relief Wrap Electrode (2" x 2")

HCPCS:

- The Relief Wrap Electrodes are made with the highest clinical standards and feature a Velcro coin so they can be easily attached to a Relief Wrap conductive brace.
- With premium white foam and a full-coverage wire connector, they are durable and safe for any user.
- Self-adhesive pad for ease of use
- Flexible carbon layer allows for maximum charge delivery and comfort
- Durable and comfortable premium white cloth backing
- Durable and reliable full-coverage wire connector
- Size: 2" x 2"
- Sealed in a tamper proof poly bag
- Velcro coin allows for easy attachment to Relief Wrap brace
- Can be reused
- For effective and high-quality electrotherapy with your Relief Wrap brace, purchase the Relief Wrap Electrodes.

PRODUCTS

Description

EP2020RW2	Electrode - Relief Wrap - 2x2 1EA/
-----------	------------------------------------