

Standard 45" TENS Lead Wire, 2 Pack

HCPCS:



- These Standard 45" TENS Lead Wires are universal replacement lead wires, compatible with most TENS unit and EMS muscle stimulators.
- TENS devices are commonly used for pain management rehabilitation, migraine pain or massaging sore muscles.
- 45" long TENS lead wires, compatible with most TENS unit and EMS muscle stimulators
- Approved for use with TENS 3000, TENS 7000, 7000 To Go, TopTENS, and InTENSity models including Twin Stim III, Twin Stim Plus 3rd Edition, Twin Stim IV, IF Combo, IF Combo II, Select Combo I and II, Micro Combo 5000 Hybrid, and InTENSity at Home
- Lead wires should be replaced every 4-6 months
- Quantity: 2 per package
- When your TENS unit's lead wire needs replacing, purchase this pack of Standard 45" Lead Wires.

PRODUCTS

Description

WW3005

Lead Wire - Standard 45 1EA/