

TheraSonic Conductive Gel, 16-oz. Bottle

HCPCS:

- The TheraSonic Ultrasound Conductive Gel is acoustically correct for a broad range of ultrasound frequencies, making it compatible with any ultrasound device.
- It is recommended for all procedures in which a viscous gel is required.
- Acoustically correct for the broad range of ultrasound frequencies
- Can be used for a variety of purposes, including fetal doppler baby heart monitor, doppler ultrasound, ultrasonic cavitation gel, and therapeutic ultrasound procedures
- Spreads evenly and won't stain clothing or damage transducers
- Water-soluble, hypoallergenic, bacteriostatic, and non-irritating
- Recommended for all procedures in which a viscous gel is required
- 16 oz.
- bottle; available in a variety of other sizes
- For an effective viscous gel for ultrasound procedures, purchase the TheraSonic Ultrasound Conductive Gel.

PRODUCTS

Description

LS5266	TheraSonic Conductive Gel, 16- 12EA/CV
--------	--
