

Support Kit

HCPCS: E1399-1,E1399-1



- Support Kit includes Trunk Harness, Abductor and Lateral Supports
- Trunk Harness:
 - Adjustable zippered trunk harness provides focused support on the sternum and rib cage
 - It's semi-stretch material allows controlled movement, prevents forward slumping and also maintains correct midline seated position
 - Straps with cam lock buckles secure the trunk harness
- Abductor:
 - Depth adjustable
 - Positions user's knee controls spastic muscles that might cause their legs to cross
- Lateral Supports:
 - Assists proper stabilization and alignment of the upper body
 - Height and width adjustable
 - Sold as a pair and each lateral is individually adjustable

PRODUCTS

Description

FC 8025 LN	Support Kit, Lg, FC4000N, 1set/bx
FC 8025N	Support Kit, Sm, FC2000N, 1set/bx

SPECIFICATIONS

Description

CA Prop 65	[WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov .]
------------	---