

## P Pod

HCPCS:



FIGURE A



FIGURE B



FIGURE C

- Provides midline alignment of trunk, pelvis, and lower extremities Simple/Lightweight design is easy for parents and caregivers to use and clean Adjustable angle based on alignment of Pod on bean bag base Provides clinically appropriate postural support through use of the molded foam liner Pelvic belt keeps user in desired position and headrest provides additional head positioning Available in 4 Sizes (Small/Medium/Large/ Extra-Large) and now available in wipe-clean Dartex fabric

### PRODUCTS

### Description

PRODUCTS	Description
APP-2000S-DX	P Pod Support, S, Wipeable Dartex
APP-2000S-MP	P Pod, S, Mermaid Pink
APP-2000S-PB	P Pod, S, Pool Blue
APP-3000M-DX	P Pod Support, M, Wipeable Dartex
APP-3000M-MP	P Pod Support, M, Mermaid Pink
APP-3000M-PB	P Pod, M, Pool Blue
APP-4000L-DX	P Pod Support, L, Wipeable Dartex
APP-4000L-JP	P Pod, L, Jellyfish Purple
APP-4000L-MP	P Pod, L, Mermaid Pink
APP-4000L-WB	P Pod, L, Whale Blue
APP-5000XL-DX	P Pod Support, XL, Wipeable Dartex
APP-5000XL-JP	P Pod, XL, Jellyfish Purple
APP-5000XL-MP	P Pod, XL, Mermaid Pink
APP-5000XL-WB	P Pod, XL, Whale Blue