

## Kanga Adult Folding Tilt-in-Space Wheelchair

HCPCS:



FIGURE A



FIGURE B



FIGURE C

- Folding frame tilt-in-space wheelchair with 45° of tilt
- Adjustable frame width and depth without additional components
- Adjustable height back canes
- Adjustable angle backrest
- Height adjustable & removable T style armrests
- Width adjustable solid seat pan included
- 20" rear wheels included 12" rear wheels optional
- Removable anti tippers included
- Transit brackets included for securing unoccupied wheelchair in vehicle
- Ideal solution for rental or long term care markets
- Must order either footrests or elevating legrests separately
- Must order seat cushion, backrest, and headrest for complete wheelchair system

### PRODUCTS

#### Description

KGA-8016-2GB	Kanga Tilt in Space, 16" W, Adj Depth
KGA-8018-2GB	Kanga Tilt in Space, 18" W, Adj Depth
KGA-8020-2GB	Kanga Tilt in Space, 20" W, Adj Depth
KGP-8010-2GR	Kanga Tilt in Space, 10" W, Adj Depth
KGP-8012-2GR	Kanga Tilt in Space, 12" W, Adj Depth
KGP-8014-2GR	Kanga Tilt in Space, 14" W, Adj Depth

### SPECIFICATIONS

#### Description

Back Angle Range	90° - 110°
Footrest Adjustment Range	11.5" - 17"
Overall Width	24" 26" 28" 12"
Tilt Angle Range	0-45°

### ACCESSORIES

#### Description

KGA-ELR-2G	Kanga Adult Elevating Legrests, Pr., EA
KGA-FR-AA	Kanga Adult Footrests, Adj Angle, Pr., EA
KG-RW12	Kanga 12" Rear Wheel Assembly, Pr., EA
NXCABF1618	Biofit Seat Cushion, 16 x 18, EA
NXCABF1818	Biofit Seat Cushion, 18 x 18, EA
NPCB1620	Classic Fixed HT DC Backrest, 16" W, EA
NPCB1820	Classic Fixed HT DC Backrest, 18" W, EA
NXCABF2020	Biofit Seat Cushion, 20 x 20, EA

**ACCESSORIES****Description**

NPCB2020	Classic Fixed HT DC Backrest, 20" W, EA
NIB-0784	Multi Fit Headrest Pad, EA
NIB-0747	Multi Axis Headrest Hardware, EA