

5" Wheel kit for WKHD-600 Walker

HCPCS:

- The 5" Wheel Kit allows your walker to glide smoothly across rough surfaces.
- They are designed for use with the WKHD-600 Bariatric Walker but may be compatible with other walkers with 1" tubing.
- 5" wheels that glide smoothly and easily across floor
- 8 adjustment holes for a wide range of height adjustments
- Compatible with WKHD-600 Walker and other walkers with 1" tubing
- Sold by the pair; 2 pairs/case
- For smoother transportation with your WKHD-600 Walker, purchase the 5" Wheel Kit.

PRODUCTS	Description
WK-HDW	5 Wheel kit for WKHD-600 and 2PR/CV