

6" STEEL ROLLATOR, BLUE, KD

HCPCS:



- The ProBasics Steel Rollator is the perfect blend of style, comfort and affordability.
- Constructed of durable powder-coated steel, it offers great maneuverability and stability.
- It features easy-to-operate ergonomic hand brakes, a removable padded backrest and a water-resistant nylon pouch.
- Equipped with an easy folding mechanism for storage and transport.
- Supports patient weight up to 300 pounds.
- The new knock-down version features a quick tool-free first time assembly.
- Durable powder-coated steel frame
- Removable curved padded backrest
- Comfortable, flip-up padded seat
- Easy-to-operate ergonomic hand brakes
- Height-adjustable for a custom fit

PRODUCTS

Description

RLS6BL-KD

6" STEEL ROLLATOR, BLUE, KD 1EA/