

Airgo MiniQuad Ultra-Stable Cane Tip

HCPCS:



FIGURE A




FIGURE B



FIGURE C



FIGURE D



FIGURE E

- The new & innovative Airgo MiniQuad ultra stable cane tip is 54% lighter and 80% more compact than a large base quad cane
- It's unique design re-centers itself for maximum stability and allows your cane to stand independently
- Fits most aluminum canes with a tubing diameter of 3/4 (1.9 cm)
- Does not fit folding canes

PRODUCTS

Description

780-030	AIRGO MINIQ, BLACK, 1EA,6CV,24CT, RET
---------	---------------------------------------

SPECIFICATIONS

Description

Warranty	Limited 90 Days
Weight Capacity	300 lbs (136 Kg)