

Bariatric Aluminum Folding Cane, Height Adjustable

HCPCS:



FIGURE A



FIGURE B



FIGURE C

- Cane folds into four convenient parts for an easy storage
- Attractive wood handle with brass collar
- Manufactured with sturdy, extruded aluminum tubing
- Comes with plastic clip to hang cane when folded
- Handle height adjusts from 33" to 37"
- Weight Capacity 500 lb

PRODUCTS

Description

| | |
|------------|---------------------------------|
| RTL10304HD | Cane,HD,Folding,Black,1EA,12/CS |
|------------|---------------------------------|

SPECIFICATIONS

Description

| | |
|-----------------|------------------|
| Warranty | Limited Lifetime |
| Handle Height | 33" - 37" |
| Weight Capacity | 500 lb |