

## Bariatric Aluminum Folding Cane, Height Adjustable

HCPCS:



FIGURE A



FIGURE B



FIGURE C

- Cane folds into four convenient parts for an easy storage
- Attractive wood handle with brass collar
- Manufactured with sturdy, extruded aluminum tubing
- Comes with plastic clip to hang cane when folded
- Handle height adjusts from 33" to 37"
- Weight Capacity 500 lb

### PRODUCTS

#### Description

|            |                                 |
|------------|---------------------------------|
| RTL10304HD | Cane,HD,Folding,Black,1EA,12/CS |
|------------|---------------------------------|

### SPECIFICATIONS

#### Description

|                 |                  |
|-----------------|------------------|
| Warranty        | Limited Lifetime |
| Handle Height   | 33" - 37"        |
| Weight Capacity | 500 lb           |