

Hugo Ergonomic Derby Handle Cane, with Claw Cane Tip

HCPCS:



FIGURE A

- Our unique Claw cane tip has a design that re-centers itself for maximum stability
- It is 57% lighter than a QuadPod, and allows your cane to stand on its own
- Claw tip is compatible with Ice-O-Grip accessory (sold separately)
- For people who could use a hand with balance while walking, our Hugo derby handle cane offers the perfect amount of stability and support
- The 12 push-button height settings ensure this cane can be adjusted to fit you properly
- The ergonomic handle with a cushioned grip helps absorb shocks
- It is ideal for those with arthritis
- Built for longevity from ultra-durable lightweight anodized aluminum
- Includes a cane strap with reflective patch for added night time safety
- Fits most people 56'5

PRODUCTS

Description

731-459	Hugo Claw Cane, Rosebud, 1EA, 6CV, Ret.
731-469	Hugo Claw Cane, Ebony, 1EA, 6CV, Ret.

SPECIFICATIONS

Description

Weight Capacity	300 lb
-----------------	--------