

ProBasics Bariatric Quad Cane, Large Base (Black), 2/ca

HCPCS:

- ProBasics Quad Canes feature a unique four-point design which allows for greater weight-bearing and stability than a single-point cane.
- The handle rotates allowing for a proper fit with right or left hand use.
- Small and large base models available in black finish with a limited lifetime warranty.
- Unique four-point design allows for greater weight bearing and stability than a single-point cane
- Non-skid rubber tips on all four feet enhance durability for long term use
- Small base design allows for easier maneuverability while the large base design provides more stability
- Handle rotates for left or right hand use
- Quad canes stand up on their own when not in use
- Easy push button height adjustment in 1-inch increments with a locking, anti-rattle ring allowing for quick, safe adjustment
- Comfortable, contoured vinyl handle with sculpted grip reduces chance of slippage
- Lightweight aluminum upper shaft with welded steel base

PRODUCTS	Description
CNQHLBB	HEAVY DUTY QUAD CANE, LG BASE 2EA/CV