

ProBasics Bariatric Two-Button Release Folding Walker, Aluminum, No Wheels,

HCPCS:

- The ProBasics Bariatric Two-Button Release Folding Walkers provide maximum strength and a rugged design, all backed by a limited lifetime warranty.
- Ample height adjustment options coupled with a two-button folding mechanism make ProBasics Folding Walkers easy to use and accommodating to the most demanding patient needs.
- Available with or without 5" wheels and in junior or adult sizes.
- Support patient weights up to 500 pounds.
- Attractive, lightweight aluminum frame
- Adjustable height settings
- Easy-to-use push button folding mechanism provides audible click when locked
- Front cross brace near top of walker allows for full stride
- Each side operates independently to allow easy movement through narrow spaces and greater stability while standing
- Vinyl contoured hand grip for added comfort and long lasting wear
- Rear glide cap allows for walker to slide easily and smoothly over most surfaces (come standard on models with 5" wheels)

PRODUCTS

Description

WKABN2B

ALUM BARIATRIC WALKER 2 BUTTON 2EA/CV